

CellVision: Microscopy Auto-Analyst with AI Narration

Buildathon Project Guide

Problem Statement

The Challenge

Cell biology researchers spend **5-10 hours per experiment** manually analyzing microscopy images - counting cells, measuring morphology, identifying phenotypes, and writing detailed figure legends for publications. This bottleneck slows research progress and introduces human error and inconsistency.

The Impact

- PhD students waste **30% of their time** on repetitive image analysis
- Manual cell counting varies **20-40% between observers**
- Writing publication-quality figure legends takes **30-60 minutes per image**
- Small labs without image analysis expertise cannot compete with well-funded institutions

Our Solution

CellVision democratizes expert-level microscopy analysis by combining computer vision with natural language AI. Upload any cell image and receive instant segmentation, quantitative metrics, and publication-ready figure legends generated by GPT-4 Vision - transforming a 2-hour manual task into a 30-second automated analysis.

Technical Roadmap (5 Hours)

Hour 0-1: Environment Setup & Model Testing

Installation Requirements

```
bash

# Core installations
pip install cellpose numpy opencv-python streamlit
pip install openai pillow matplotlib scikit-image
```

Initial Model Test

```
python
```

```
# Download and test CellPose
from cellpose import models
import numpy as np
from skimage import io

# Initialize model (use CPU for reliability)
model = models.Cellpose(gpu=False, model_type='cyto2')

# Test on sample image
test_img = io.imread('sample_cells.png')
masks, flows, styles, diams = model.eval(test_img)
print(f"Detected {len(np.unique(masks))-1} cells") # Verify it works
```



Deliverable: Working CellPose installation with successful test segmentation

Hour 1-2: Core Analysis Pipeline

Main Analysis Function

```
python
```

```

from skimage.measure import regionprops_table
import numpy as np

def analyze_microscopy_image(image_path):
    """
    Complete microscopy image analysis pipeline.

    Args:
        image_path: Path to microscopy image

    Returns:
        masks: Segmentation masks for each cell
        metrics: Dictionary of quantitative measurements
    """
    # 1. Load and preprocess
    img = io.imread(image_path)

    # 2. Cell segmentation with CellPose
    masks, _, _, diams = model.eval(img, diameter=30)

    # 3. Extract quantitative metrics
    cell_count = len(np.unique(masks)) - 1

    # 4. Morphology analysis
    props = regionprops_table(masks, img,
                             properties=['area', 'perimeter', 'eccentricity', 'solidity'])

    # 5. Calculate statistics
    metrics = {
        'total_cells': cell_count,
        'avg_area': np.mean(props['area']) if len(props['area']) > 0 else 0,
        'avg_circularity': np.mean(4*np.pi*props['area']/props['perimeter']**2) if len(props['area']) > 0 else 0,
        'density': cell_count / (img.shape[0] * img.shape[1]),
        'size_variation': np.std(props['area']) / np.mean(props['area']) if len(props['area']) > 0 else 0
    }

    return masks, metrics

```

✓ **Deliverable:** Function that inputs image → outputs segmentation masks + metrics dictionary

Hour 2-3: GPT-4 Vision Integration

AI Narrative Generation

python

```

import base64
from openai import OpenAI
import matplotlib.pyplot as plt
from skimage import io

def generate_analysis_narrative(image_path, masks, metrics, api_key):
    """
    Generate publication-quality figure legend using GPT-4 Vision.

    Args:
        image_path: Path to original image
        masks: Segmentation masks from CellPose
        metrics: Quantitative metrics dictionary
        api_key: Azure OpenAI API key

    Returns:
        str: Publication-ready figure legend
    """
    # 1. Create visualization
    fig, (ax1, ax2) = plt.subplots(1, 2, figsize=(12, 6))
    ax1.imshow(io.imread(image_path))
    ax1.set_title('Original')
    ax1.axis('off')
    ax2.imshow(masks, cmap='tab20')
    ax2.set_title(f'Segmented: {metrics["total_cells"]} cells')
    ax2.axis('off')
    plt.savefig('analysis.png', dpi=100, bbox_inches='tight')
    plt.close()

    # 2. Encode image for GPT-4 Vision
    with open('analysis.png', 'rb') as f:
        img_base64 = base64.b64encode(f.read()).decode('utf-8')

    # 3. Call GPT-4 Vision with metrics context
    client = OpenAI(api_key=api_key)

    prompt = f"""You are an expert cell biologist. Analyze this microscopy image.

    Quantitative metrics detected:
    - Cell count: {metrics['total_cells']}
    - Average cell area: {metrics['avg_area']:.1f} pixels²
    - Cell density: {metrics['density']:.4f} cells/pixel²
    - Size variation coefficient: {metrics['size_variation']:.2f}
    """

```

```
- Average circularity: {metrics['avg_circularity']:.2f}
```

Generate a publication-quality figure legend describing:

1. Cell morphology and distribution patterns
2. Notable features or abnormalities
3. Quantitative summary with the provided statistics

Format as: "Figure X: [Comprehensive scientific description]"

Keep it concise but scientifically precise.

```
"""
```

```
response = client.chat.completions.create(  
    model="gpt-4-vision-preview",  
    messages=[  
        {"role": "user",  
         "content": [  
             {"type": "text", "text": prompt},  
             {"type": "image_url", "image_url": {  
                 "url": f"data:image/png;base64,{img_base64}"  
             }}  
         ]  
    },  
    ],  
    max_tokens=500  
)  
  
return response.choices[0].message.content
```

✓ **Deliverable:** Function that generates natural language analysis from image + metrics

Hour 3-4: Streamlit Interface

Complete Web Application

```
python
```

```

import streamlit as st
import matplotlib.pyplot as plt
from datetime import datetime
import os

def main():
    # Page config
    st.set_page_config(
        page_title="CellVision",
        page_icon="🔬",
        layout="wide"
    )

    # Header
    st.title("🔬 CellVision: AI-Powered Microscopy Analysis")
    st.markdown("Transform microscopy images into quantitative insights and publication-ready descriptions in seconds.")

    # Sidebar for API key
    with st.sidebar:
        st.header("⚙️ Configuration")
        api_key = st.text_input("OpenAI API Key", type="password")
        st.markdown("---")
        st.markdown("### 🖼️ Sample Images")
        if st.button("Load Cancer Cell Example"):
            # Load pre-saved example
            pass

    # Main content
    col1, col2 = st.columns([1, 1])

    with col1:
        st.subheader("📁 Upload Image")
        uploaded_file = st.file_uploader(
            "Choose a microscopy image",
            type=['png', 'jpg', 'jpeg', 'tif', 'tiff'],
            help="Supported formats: PNG, JPEG, TIFF"
        )

        if uploaded_file:
            # Save uploaded file
            temp_path = f"temp_{datetime.now().timestamp()}.png"
            with open(temp_path, "wb") as f:
                f.write(uploaded_file.getbuffer())

```

```
# Display original
```

```
st.image(temp_path, caption="Original Image", use_column_width=True)
```

```
with col2:
```

```
st.subheader("🔍 Analysis Results")
```

```
if uploaded_file and st.button("🚀 Analyze Image", type="primary"):
```

```
    if not api_key:
```

```
        st.error("Please enter your OpenAI API key in the sidebar.")
```

```
    else:
```

```
        # Progress indicator
```

```
        progress_bar = st.progress(0)
```

```
        status_text = st.empty()
```

```
        # Step 1: Segmentation
```

```
        status_text.text("Segmenting cells...")
```

```
        progress_bar.progress(33)
```

```
        masks, metrics = analyze_microscopy_image(temp_path)
```

```
        # Display segmentation
```

```
        fig, ax = plt.subplots(figsize=(8, 8))
```

```
        ax.imshow(masks, cmap='tab20')
```

```
        ax.set_title(f"Detected {metrics['total_cells']} cells")
```

```
        ax.axis('off')
```

```
        st.pyplot(fig)
```

```
        plt.close()
```

```
        # Step 2: Metrics
```

```
        status_text.text("Calculating metrics...")
```

```
        progress_bar.progress(66)
```

```
        # Display metrics in columns
```

```
        metric_cols = st.columns(4)
```

```
        metric_cols[0].metric("Cell Count", metrics['total_cells'])
```

```
        metric_cols[1].metric("Avg Area", f"{metrics['avg_area']:.0f} px2")
```

```
        metric_cols[2].metric("Density", f"{metrics['density']:.4f}")
```

```
        metric_cols[3].metric("Circularity", f"{metrics['avg_circularity']:.2f}")
```

```
        # Step 3: AI Analysis
```

```
        status_text.text("Generating AI analysis...")
```

```
        progress_bar.progress(90)
```

```
narrative = generate_analysis_narrative(temp_path, masks, metrics, api_key)
```



```

# Complete
progress_bar.progress(100)
status_text.text("Analysis complete!")

# Display narrative
st.markdown("### 📄 Publication-Ready Figure Legend")
st.info(narrative)

# Download button
st.download_button(
    label="📎 Download Analysis Report",
    data=narrative,
    file_name=f"figure_legend_{datetime.now().strftime('%Y%m%d_%H%M%S')}.txt",
    mime="text/plain"
)

# Cleanup
os.remove(temp_path)

if __name__ == "__main__":
    main()

```

✅ **Deliverable:** Complete working web interface with upload → analyze → download flow

Hour 4-5: Testing, Polish & Demo Prep

Testing Checklist

- ☐ Test with 5 different cell types from LIVECell dataset
- ☐ Verify cell counts within 10% of manual annotation
- ☐ Ensure GPT-4 narratives are scientifically accurate
- ☐ Test edge cases: dense clusters, low contrast, artifacts
- ☐ Verify error handling for corrupted images
- ☐ Test with different image sizes (512x512 to 2048x2048)

Polish Tasks

python

Add these enhancements if time permits:

```
def add_batch_processing():
    """Allow multiple image upload and analysis"""
    uploaded_files = st.file_uploader("Choose images", accept_multiple_files=True)
    if uploaded_files:
        for file in uploaded_files:
            # Process each file
            pass

def export_pdf_report():
    """Generate comprehensive PDF report with all results"""
    from reportlab.lib.pagesizes import letter
    from reportlab.pdfgen import canvas
    # Create PDF with images, metrics, and narrative
    pass

def add_comparison_mode():
    """Side-by-side comparison of two conditions"""
    col1, col2 = st.columns(2)
    # Allow comparison of treated vs untreated cells
    pass
```

Demo Data Preparation

1. Download test images:

- Healthy epithelial cells
- Cancer cells (HeLa)
- Drug-treated cells showing apoptosis

2. Create demo narrative:

- "Watch CellVision identify apoptotic cells invisible to the untrained eye"
- Show 3-step progression: healthy → cancer → treated

3. Backup plan:

- Save screenshots of each step
 - Pre-generate example outputs
 - Have PDF report ready to show
-

Team Task Assignments

Person 1 - Data Curator

- ☐ Download 20 diverse cell images from LIVECell dataset
- ☐ Organize by cell type (A172, BT474, BV2, Huh7, MCF7, etc.)
- ☐ Create metadata spreadsheet with image properties
- ☐ Select 3 best images for demo

Person 2 - Validator

- ☐ Manually count cells in 5 test images
- ☐ Create ground truth spreadsheet
- ☐ Calculate accuracy metrics
- ☐ Document edge cases where algorithm struggles

Person 3 - UX Designer

- ☐ Create Figma mockup of interface
- ☐ Design CellVision logo
- ☐ Choose color scheme (suggest: teal/purple for biotech feel)
- ☐ Create slide deck template (5 slides max)

Person 4 - Demo Narrator

- ☐ Write 90-second pitch script
 - ☐ Prepare problem → solution → impact storyline
 - ☐ Create slide deck content:
 - Slide 1: Problem (30 sec)
 - Slide 2: Solution (30 sec)
 - Slide 3: Live Demo (60 sec)
 - Slide 4: Impact & Next Steps (20 sec)
 - ☐ Practice transitions between speakers
-





Success Metrics

Technical Goals

- ☒ Process any microscopy image in **<30 seconds**
- ☒ Cell count accuracy **within 90%** of manual annotation
- ☒ Generate coherent, scientific figure legends
- ☒ Support images from **512x512 to 2048x2048** pixels

-  Handle **PNG, JPEG, and TIFF** formats

Demo Goals

-  Live demo works flawlessly with 3 different cell types
 -  Clear value proposition: **2 hours** → **30 seconds**
 -  Show quantitative improvement in consistency
 -  Generate "wow" moment with AI-generated narrative
-

Winning Pitch Framework

Opening Hook (15 seconds)

"Every day, 50,000 biology PhD students spend 2 hours counting cells by hand, clicking on dots one by one. Meanwhile, critical drug discoveries are delayed because small labs can't afford \$100K image analysis software."

Problem Validation (15 seconds)

"Manual cell counting varies 40% between researchers. One image takes 30 minutes to analyze properly. Publication-quality figure legends take another 30 minutes to write."

Solution Demo (60 seconds)

1. **Upload** cancer cell image (5 sec)
2. **Show** instant segmentation with 234 cells detected (10 sec)
3. **Display** quantitative metrics dashboard (10 sec)
4. **Reveal** AI-generated figure legend (20 sec)
5. **Compare** to manual analysis time/accuracy (15 sec)

Impact Statement (20 seconds)

"CellVision democratizes advanced microscopy analysis. Any researcher with a phone camera microscope can now produce Nature-quality image analysis. We transform a 2-hour expert task into 30 seconds of automated precision."

Ask & Next Steps (10 seconds)

"We're seeking clinical research partners to validate CellVision across 10 cell types and apply for NIH SBIR funding to build the definitive microscopy analysis platform."

Common Pitfalls to Avoid

1. **Don't oversell accuracy** - Be honest about 90% accuracy vs 100%

2. **Don't ignore edge cases** - Acknowledge when cells are too dense/overlapping
 3. **Don't skip validation** - Always compare to ground truth
 4. **Don't forget citations** - Credit CellPose and LIVECell dataset
 5. **Don't overcomplicate** - Simple, working demo beats complex, buggy features
-



Resources & References

- **CellPose Documentation:** <https://cellpose.readthedocs.io/>
 - **LIVECell Dataset:** <https://sartorius-research.github.io/LIVECell/>
 - **OpenAI Vision API:** <https://platform.openai.com/docs/guides/vision>
 - **Streamlit Docs:** <https://docs.streamlit.io/>
 - **Sample Microscopy Images:** <https://bbbc.broadinstitute.org/>
-



Why This Wins

1. **Visual Impact:** Transformation is immediately visible and impressive
2. **Real Problem:** Every biology lab faces this challenge daily
3. **Technical Innovation:** First to combine segmentation + GPT-4 Vision narration
4. **Feasible in 5 Hours:** Core functionality proven in 2 hours, polish in remaining 3
5. **Democratizes Science:** Levels playing field between large and small labs
6. **Clear Metrics:** 2 hours → 30 seconds, 40% variation → 10% variation

Remember: Judges want innovation and out-of-the-box thinking. CellVision delivers both while solving a real, painful problem that every biologist understands.